



UTAH MEDICAL PRODUCTS INC.
 7043 South 300 West • Midvale, Utah 84047
 800-533-4984 • 801-566-1200
www.utahmed.com

PRESORT STANDARD
 U. S. POSTAGE
PAID
 PERMIT NO. 3280
 SALT LAKE CITY, UT

Nutri-Lok™ ...



Utah Medical Products, Inc. has augmented Nutri-Cath[®] with the Nutri-Lok[™] adaptor.

Nutri-Lok is a locking enteral feeding-only extension set that:

- * Will not mate with a female luer IV line connector
- * Ensures a secure connection that will not accidentally slip apart
- * Is compatible via slip fit with other oral dose enteral feeding catheters, but only LOCKS with Nutri-Cath
- * Avoids inadvertent delivery of enteral feeding through the intravenous route



The oral dose Nutri-Cath hubs now have threads which allow extension set locking to the catheter hub..




Competitive sets connect via slip fit - when enteral solutions are administered, the connection becomes slippery and can easily disconnect.

why take chances?



With Nutri-Lok, a secure locking connection is intended to ensure no accidental disconnects



Nutri-Lok		Packaged 50/bx		
60" DEHP-Free Extension Set		4150007		
Nutri-Cath		Packaged 10/bx		
Oral Dose Hub	3.5 Fr	5.0 Fr	6.5 Fr	8.0 Fr
40 cm	4153517	4155017	4156517	4158017
58 cm	4153527	4155027	4156527	4158027
93 cm	4153537	4155037	4156537	4158037

For more information or to request a sample, visit:

www.utahmed.com

or call:

1-800-533-4984

UTAH MEDICAL
PRODUCTS INC.

