

# OptiSpec® Gynecology Light

Hands-free ultra-bright illumination provides visualization without inconvenience.



Utah Medical Products' patented\* OptiSpec Light is a new concept in non-colposcopic illumination of the cervix. An ultra-bright LED selected to provide a pure white light spectrum has been mounted in a small, clip-on disposable package. The result is excellent illumination of the cervix with a device that otherwise seems like it's not even there!

Cervical visualization through a colposcope with a bright white light provides critical visual information with minimal clinician fatigue. However, the use of the colposcope for other every-day exams is impractical. Other methods of cervical illumination emit a dull yellow, low intensity light, and usually require one hand to actively hold the lighting device.



- Compact light clips on to most common vaginal specula
- Unobtrusive configuration remains out of visual and working field
- Provides a simple, hands-free light that improves visualization during:
  - gyn exams
  - pap smears
  - LETZ® procedures
  - diagnosis of abnormal obstetric bleeding
  - ER exams for vaginal trauma
  - any other directly visualized vaginal procedures
- Efficient light-emitting diode (LED) provides truer color visualization with a light that is whiter than halogen bulbs

## Order Information

- Item Number: LITE-WS
- Packaging: 25 per box  
Single Use, Sterile

## Specifications

- Housing: ABS
- Light source: Ultrabright white LED

The OptiSpec Light is LATEX FREE

\*U.S. Patent 7,631,981



## United States

Utah Medical Products, Inc.  
7043 South 300 West  
Midvale, Utah 84047  
800.533.4984 ♦ 801.566.1200  
fax 801.566.2062

## Europe

Utah Medical Products Ltd.  
Athlone Business & Technology Park  
Dublin Road  
Athlone, County Westmeath, Ireland  
353.90.647.3932

[www.utahmed.com](http://www.utahmed.com)

pn 58255 rev 111111