



SECURE
~~SLIPPERY~~
WHEN WET...

DEHP
Free

Nutri-Lok™

(patent pending)

Even if you already meet JCAHO requirements for all oral enteral feeding systems, are you still taking unnecessary chances? Competitive oral sets and syringes connect via slip fit - when enteral solutions are administered, the connections may become slippery and can easily disconnect. With Nutri-Lok locking extension sets and syringes, the secure locking connections are designed to reduce accidental disconnects.

Nutri-Lok is a family of locking enteral feeding-only extension sets and syringes that:

- * Will not mate with a female luer IV line connector
- * Ensures a secure connection that will not accidentally slip apart
- * Are compatible via slip fit with other oral dose enteral feeding catheters, but only LOCK with Nutri-Cath®

unnecessary
why take chances?

Nutri-Lok

Packaged 50/bx

4150007	60" Locking Extension Set, Female Standard Luer to Male Oral
4150017	60" Locking Extension Set, All-Oral Female to Male
41560	60mL Sterile Nutri-Lok Oral Dose Syringe, individually packaged

Nutri-Cath

Packaged 10/bx

Oral Dose Hub	3.5 Fr	5.0 Fr	6.5 Fr	8.0 Fr
40 cm	4153517	4155017	4156517	4158017
58 cm	4153527	4155027	4156527	4158027
93 cm	4153537	4155037	4156537	4158037

United States

Utah Medical Products, Inc.
7043 South 300 West
Midvale, Utah 84047
Toll-free 800 533-4984
Tel 801 566-1200
Fax 801 566-2062
www.utahmed.com

UTAH MEDICAL
PRODUCTS INC.



Europe

Utah Medical Products, Ltd.
Athlone Business & Technology Park
Athlone, County Westmeath, Ireland
+353-90-647-3932