


Uri-Cath™

Closed, Pre-assembled
Urinary Drainage System

Because babies in the NICU *require*
gentle urinary drainage

Accurate measurement of urinary output
for neonates and pediatric patients

Soft silicone material minimizes tissue
irritation and encrustation



- Pre-assembled, sterile system reduces risk of infection and reduces valuable nursing time
- Stopcock with sampling access site allows for non-contaminated midstream specimen sampling
- 150 ml burette for accurate urine output measurement

Utah Medical Products, Inc
manufactures a broad range of
neonatal critical care devices
including the Gesco™ family
of products.



UTAH MEDICAL
PRODUCTS INC.





Uri-Cath™

Product Description

The Uri-Cath™ tray consists of a pre-connected silicone urinary drainage catheter, three-way stopcock for sampling purposes, and a 150 ml burette for accurate urine output measurement.

Design Features

- Biocompatible silicone tubing
- Encased radiopaque stripe
- Soft, rounded tip
- Stopcock with latex-free injection site
- 150 ml collection burette with drainage stopcock
- 3 distal openings
- Available in three sizes:
3.5 French, 5.0 French, 8.0 French

Ordering Information

Urinary Drainage Trays

Model #	Description	Unit of Measure	Image
4193507	3.5 French Uri-Cath™ set with silicone catheter	10/BOX	
4195007	5.0 French Uri-Cath™ set with silicone catheter	10/BOX	
4198007	8.0 French Uri-Cath™ set with silicone catheter	10/BOX	

Urinary Drainage Catheters

Model #	Description	Unit of Measure	Image
4193505	3.5 French Uri-Cath™ catheter only	10/BOX	
4195005	5.0 French Uri-Cath™ catheter only	10/BOX	
4198005	8.0 French Uri-Cath™ catheter only	10/BOX	

UNITED STATES
7043 South 300 West
Midvale, Utah 84047

(801) 566-1200
Toll Free 1-800-533-4984
Fax (801) 566-2062

www.utahmed.com

EUROPE EC REP
Utah Medical Products, Ltd.
Athlone Business & Technology Park
Athlone, County Westmeath, Ireland
+353-90-647-3932



Gesco™

P/N 58276 rev 031108