



# NUTRI-LOK<sup>®</sup>


## ENTERAL FEED

### 500ml Pump Set

~with Nutri-Lok Connector

DEHP  
Free

- Oral dose locking connection ensures bag and Nutri-Cath catheter will not accidentally slip apart
- Incorporates an Anti-Free Flow Valve as a safety feature to prevent free flow of formula into the patient
- Nutri-Lok connector will not mate with a female luer IV connector
- DEHP-Free, Latex-Free



Nutri-Cath<sup>®</sup> silicone  
feeding catheter



For use with KANGAROO ePump<sup>™</sup> enteral feeding pumps

UTAH MEDICAL  
PRODUCTS INC.



VISIT [www.utahmed.com/nutricath.htm](http://www.utahmed.com/nutricath.htm) OR CALL TOLL FREE  
800-533-4984 FOR CUSTOMER SERVICE

When you're ready to become JCAHO compliant<sup>1</sup>,

UTAH MEDICAL PRODUCTS INC.



Utah Medical Products, Inc. provides a complete enteral-only system for your patient safety needs.

## 1- Syringe feedings

- via gravity feed:

- \*Secure locking connections prevent inadvertent disconnections from occurring.
- \*Clinician's valuable time will not be wasted cleaning up bedside spills.
- \*Mother's precious breast milk will be entirely delivered to baby.



- via syringe pump:

- \*Monoject™ oral dispensing syringes are recognized by existing syringe pumps, eliminating potential dosage errors.
- \*Multiple options for Nutri-Lok extension sets allow hospitals to choose a continuous feed adaptor that fits with each unit's unique feeding protocol.



## 2- Bag feedings

- \*Nutri-Lok Enteral Feed Bags provide patients who require larger feeds with the clinical benefits and safety advantages of the Nutri-Lok feeding family.

<sup>1</sup> [http://www.jointcommission.org/SentinelEvents/SentinelEventAlert/sea\\_36.htm](http://www.jointcommission.org/SentinelEvents/SentinelEventAlert/sea_36.htm)

	Nutri-Cath feeding catheter with Nutri-Lok				Packaged 10/bx
	3.5 Fr	5.0 Fr	6.5 Fr	8.0 Fr	
40 cm	4153517	4155017	4156517	4158017	
58 cm	4153527	4155027	4156527	4158027	
93 cm	4153537	4155037	4156537	4158037	

Nutri-Lok Continuous Feed Adaptors		Packaged 50/bx
4150017	60" Locking Ext Set, Female Nutri-Lok to Male Nutri-Lok	
4150027	60" Locking Ext Set w/ side port, Female Nutri-Lok to Male Nutri-Lok	
4150037	60" Universal Ext Set w/ side port, Female Nutri-Lok to Stepped Adaptor	

Enteral Feed Bag		Packaged 30/bx
4150047	Nutri-Lok Enteral Feed 500ml Pump Set with Nutri-Lok Connector for use with Kangaroo ePump™ feeding pumps	

Sterile Nutri-Lok oral dispensing syringe, individually packaged		
Volume	Part #	
1 mL:	41501	100/bx
3 mL:	41503	100/bx
6 mL:	41506	100/bx
12 mL:	41512	100/bx
20 mL:	41520	100/bx
35 mL:	41535	100/bx
60 mL:	41560	50/bx

Sterile Nutri-Lok oral dispensing syringe caps, 5 per package		
Caps	Part #	
	415CAP	200/bx

Sterile Nutri-Lok oral dose straws, individually packaged		
Straws	Part #	
	415STR	50/bx

**United States**  
 Utah Medical Products, Inc.  
 7043 South 300 West  
 Midvale, Utah 84047  
 Toll-free 800 533-4984  
 Tel 801 566-1200  
 Fax 801 566-2062

**Europe**  
 Utah Medical Products, Ltd.  
 Athlone Business & Technology Park  
 Athlone, County Westmeath, Ireland  
 +353-90-647-3932

[www.utahmed.com](http://www.utahmed.com)