

LETZ[®] Coated Instruments

Quality, Durability, and Selection

Utah Medical Products, Inc. offers LETZ-coated instruments that are the ideal combination of quality stainless steel vaginal specula and extremely durable nylon coating.

A Broad Range of Sterilization Options


These LETZ instruments have been qualified to withstand sterilization by:

- Autoclave (up to 280°F)
- Ethylene Oxide
- Sterrad*

A Broad Range of Configurations

Patient anatomy and provider preference dictate the need for different product configurations. Utah Medical Products provides choices in blade style (Graves, Pederson, etc.) and length (from small to extra large), access (introital openings ranging from standard to extra wide widths and open-sided), and even an option with built-in lateral wall retractors for collapsing walls.

A Comparison of Specula



*Sterrad is a registered trademark of Advanced Sterilization Products, Division of Ethicon Inc, a Johnson and Johnson company

- Integration of advanced materials and coating process produces an extremely durable coating that stays intact for long instrument life.
- Built-in smoke tube on Graves and Pederson models eliminate the need for "DSE" tubing, reducing procedure cost.
- Certified to be compatible with Sterrad process.



UTAH MEDICAL PRODUCTS INC.

United States

Utah Medical Products, Inc.
7043 South 300 West
Midvale, Utah 84047
800.533.4984 • 801.566.1200
fax 801.566.2062

Ireland

Utah Medical Products Ltd.
Athlone Business & Technology Park
Athlone, County Westmeath
(090) 647.3932 • fax (090) 647.5608

United Kingdom

Femcare-Nikomed Ltd.
Stuart Court, Spursholt
Place, Salisbury Road
Romsey, Hampshire SO51 6DJ
(01794) 525.100 • fax (01794) 525.101

www.utahmed.com

LETZ® Coated Instruments

Blade lengths for Graves specula, Pederson specula, and Four-Way Expander

- Small 3½" 9cm
- Medium 4¼" 11cm
- Large 4¾" 12cm
- Extra Long 6" 15cm

Order Information

Graves Specula

Standard "duckbill" rounded blade shape for normal use

ESI-100	Standard, Small
ESI-101	Standard, Medium
ESI-102	Standard, Large
ES-15162-XSLT	Standard, Extra Long
ESI-151	Wide, Medium
ESI-152	Wide, Large
ESI-171	View-Maxi Speculum, Medium
ESI-172	View-Maxi Speculum, Large
ESI-131	Weisman-Graves, Left Side Open, Medium
ESI-132	Weisman-Graves, Left Side Open, Large
ESI-133	Weisman-Graves, Right Side Open, Medium
ESI-134	Weisman-Graves, Right Side Open, Large



Pederson Specula

Straight, narrow blade shape for narrow introitus

ESI-110	Standard, Small
ESI-111	Standard, Medium
ESI-112	Standard, Large
ESI-117	View-Maxi, Medium
ESI-118	View-Maxi, Large



Four-Way Expanders

Manage collapsing vaginal walls with a single instrument

ES-16122-MLE	Four Way Vaginal Expander, Medium
ES-16132-LLE	Four Way Vaginal Expander, Large
ES-16135-XLE	Four Way Vaginal Expander, Extra Large

also available in uncoated stainless steel and autoclavable resin models



ES-16132-LLE
Large LETZ-Coated
Four-Way Expander



United States

Utah Medical Products, Inc.
7043 South 300 West
Midvale, Utah 84047
800.533.4984 • 801.566.1200
fax 801.566.2062

Ireland

Utah Medical Products Ltd.
Athlone Business & Technology Park
Athlone, County Westmeath
(090) 647.3932 • fax (090) 647.5608

United Kingdom

Femcare-Nikomed Ltd.
Stuart Court, Spursholt
Place, Salisbury Road
Romsey, Hampshire SO51 6DJ
(01794) 525.100 • fax (01794) 525.101

www.utahmed.com