


# DISPOSABLE-A-HOOD™

## INFANT OXYGEN HOOD



*Designed to promote precise control of Neutral Thermal Environment*

*Ready-to-use design saves valuable time*

*Disposability eliminates risk of cross-contamination*

*Inlet port diffuser prevents oxygen from flowing directly onto baby*

*New extra large size*

*Convenient top port accommodates temperature and oxygen sensors*


*Unique "laminar flow" design maintains consistent FiO<sub>2</sub> levels and improves CO<sub>2</sub> flushing*

*Padded neck opening is gentle to baby's skin*

*Raised base ports allow for easy placement of tubing and cables*

*Clear, seamless material allows total visibility of neonate from all angles*

*Weighted, reusable collar ensures hood stability and promotes ease of use*



The Disposa-Hood™ Disposable Infant Oxygen Hood offers the most popular advantages found in reusable hoods. Distortion-free neonate visualization and excellent FiO<sub>2</sub> distribution are integrated in a disposable hood, so cross contamination is no longer a risk. With Disposa-Hood, the toils of inter-use cleaning, physical tracking, and regular bioburden buildup testing are eliminated, conserving the hospital's valuable time and resources. Also eliminated are the unsightly scratches that accumulate with repeated use and cleaning, as well as the expensive replacement cost for a damaged reusable hood.

Disposa-Hood also provides better performance over cumbersome, inefficient tents. The Disposa-Hood is ready to use right out of the box, saving the clinician the time and frustration involved with setting up the tent. The seamless, transparent design eliminates guesswork when the neonate requires visual assessment. Efficient distribution of the FiO<sub>2</sub> and flushing of CO<sub>2</sub> are important benefits of the rounded Disposa-Hood shape, a flow design which cannot be found in the available tents.

With the wide range of sizes available and the advantages that Disposa-Hood offers, there is no reason to accept less than the best possible oxygen delivery system.

## DISPOSA-HOOD INFORMATION

<i>Catalog Number</i>	<i>Size</i>	<i>Overall Dimensions</i>			<i>Neck Opening</i>		<i>Quantity</i>
		<i>Length</i>	<i>Width</i>	<i>Height</i>	<i>Width</i>	<i>Height</i>	
0300	Medium	7.0"	7.0"	5.4"	3.9"	2.5"	24/case
0301	Small	7.0"	7.0"	5.4"	2.5"	1.8"	24/case
0305	Large	8.0"	8.3"	6.9"	5.0"	3.9"	16/case
0305W	Large, Wide Neck	8.0"	8.3"	6.9"	5.0"	2.9"	16/case
0313	Extra Large	13.0"	10.4"	8.5"	6.3"	5.0"	12/case

## STARTER KIT

<i>Catalog Number</i>	<i>Description</i>	<i>Quantity</i>
0325	Disposa-Hood Starter Kit Contains one each 0300, 0301, 0305, and 0305W hoods, plus one each 0302 and 0308 collars	6 pcs

## REUSABLE COLLARS

<i>Catalog Number</i>	<i>Description</i>	<i>Quantity</i>
0302	Collar for Small/Medium Disposa-Hoods	2/box
0308	Collar for Large Disposa-Hoods	2/box
0314	Collar for Extra Large Disposa-Hood	2/box

UTAH MEDICAL  
PRODUCTS INC.



### United States

7043 South 300 West  
Midvale, Utah 84047  
Tel: 801/566-1200  
Toll Free: 800/533-4984  
Fax: 801/566-2062

### Europe

Utah Medical Products, Ltd.  
Athlone Business & Technology Park  
Athlone, County Westmeath, Ireland  
+ 353-90-647-3932

READY TO USE  
EFFICIENT  
TRANSPARENT  
DISPOSABLE  
CONVENIENT  
PORTABLE  
LAMINAR FLOW